

Inhibitor of nuclear factor kappa-B kinase subunit beta

Catalog No: #AP84393



Package Size: #AP84393-1 50ug #AP84393-2 100ug #AP84393-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Inhibitor of nuclear factor kappa-B kinase subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	I-kappa-B-kinase beta,I-kappa-B kinase 2,IKK2,Nuclear factor NF-kappa-B inhibitor kinase beta,NFKB1KB,Ikkb
Accession No.	Uniprot ID:Q9QY78
Uniprot	Q9QY78
GeneID	84351;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ikkb

Note: This product is for in vitro research use only