

# Interleukin-2 receptor subunit alpha

Catalog No: #AP84395



Package Size: #AP84395-1 50ug #AP84395-2 100ug #AP84395-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Interleukin-2 receptor subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	IL-2 receptor subunit alpha,CD25
Accession No.	Uniprot ID:O02733
Uniprot	O02733
GeneID	396814;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:IL2RA

Note: This product is for in vitro research use only