

## Dyslexia-associated protein KIAA0319-like protein

Catalog No: #AP84398



Package Size: #AP84398-1 50ug #AP84398-2 100ug #AP84398-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Dyslexia-associated protein KIAA0319-like protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Adeno-associated virus receptor,AAVR,AAVR,KIAA1837,PP791
Accession No.	Uniprot ID:Q8IZA0
Uniprot	Q8IZA0
GeneID	79932;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: KIAA0319L

Note: This product is for in vitro research use only