

High mobility group protein B2

Catalog No: #AP84401



Package Size: #AP84401-1 50ug #AP84401-2 100ug #AP84401-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	High mobility group protein B2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	High mobility group protein 2,HMG-2,Hmg2
Accession No.	Uniprot ID:P30681
Uniprot	P30681
GeneID	97165;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hmgb2

Note: This product is for in vitro research use only