

Histamine H4 receptor

Catalog No: #AP84404



Package Size: #AP84404-1 50ug #AP84404-2 100ug #AP84404-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Histamine H4 receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	H4R
Accession No.	Uniprot ID:Q91ZY2
Uniprot	Q91ZY2
GeneID	225192;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Hrh4

Note: This product is for in vitro research use only