

Serine protease inhibitor Kazal-type 1

Catalog No: #AP84405



Package Size: #AP84405-1 50ug #AP84405-2 100ug #AP84405-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Serine protease inhibitor Kazal-type 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	Pancreatic secretory trypsin inhibitor,PSTI
Accession No.	Uniprot ID:P00998
Uniprot	P00998
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:SPINK1

Note: This product is for in vitro research use only