

Integrin beta-1

Catalog No: #AP84414



Package Size: #AP84414-1 50ug #AP84414-2 100ug #AP84414-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Integrin beta-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Beta oligodendroglia,Beta OL,Fibronectin receptor subunit beta,VLA-4 subunit beta,CD29
Accession No.	Uniprot ID:P49134
Uniprot	P49134
GeneID	24511;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Itgb1

Note: This product is for in vitro research use only