

Protein kinase C alpha type

Catalog No: #AP84420



Package Size: #AP84420-1 50ug #AP84420-2 100ug #AP84420-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein kinase C alpha type
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PKC-A,PKCA,PRKACA
Accession No.	Uniprot ID:P17252
Uniprot	P17252
GeneID	5578;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PRKCA

Note: This product is for in vitro research use only