

# Glutamate receptor 1

Catalog No: #AP84423

Package Size: #AP84423-1 50ug #AP84423-2 100ug #AP84423-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Glutamate receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GluR-1,AMPA-selective glutamate receptor 1,GluR-A,GluR-K1,Glutamate receptor ionotropic, AMPA 1,GluA1,GLUH1,GLUR1
Accession No.	Uniprot ID:P42261
Uniprot	P42261
GeneID	2890;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:GRIA1

Note: This product is for in vitro research use only