

## Inter alpha-trypsin inhibitor, heavy chain 4

Catalog No: #AP84425



Package Size: #AP84425-1 50ug #AP84425-2 100ug #AP84425-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Inter alpha-trypsin inhibitor, heavy chain 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	ITI heavy chain H4
Accession No.	Uniprot ID:A6X935
Uniprot	A6X935
GeneID	16427;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Itih4

Note: This product is for in vitro research use only