

Potassium voltage-gated channel subfamily B member 2



Catalog No: #AP84427

Orders: order@signalwayantibody.com

Package Size: #AP84427-1 50ug #AP84427-2 100ug #AP84427-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Potassium voltage-gated channel subfamily B member 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CDRK,Voltage-gated potassium channel subunit Kv2.2
Accession No.	Uniprot ID:Q63099
Uniprot	Q63099
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Kcnb2

Note: This product is for in vitro research use only