

Galectin-3-binding protein

Catalog No: #AP84429



Package Size: #AP84429-1 50ug #AP84429-2 100ug #AP84429-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Galectin-3-binding protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cyp-C-associated protein,CyCAP,Lectin galactoside-binding soluble 3-binding protein,Mac-2-binding protein,MAC2BP,Protein 90K
Accession No.	Uniprot ID:O70513
Uniprot	O70513
GeneID	245955;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lgals3bp

Note: This product is for in vitro research use only