

Interleukin-24

Catalog No: #AP84430



Package Size: #AP84430-1 50ug #AP84430-2 100ug #AP84430-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-24
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IL-24,IL-4-induced secreted protein,Melanoma differentiation-associated gene 7 protein,MDA-7,Th2-specific cytokine FISP,Mda7
Accession No.	Uniprot ID:Q925S4
Uniprot	Q925S4
GeneID	93672;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL24

Note: This product is for in vitro research use only