

Interleukin-2 receptor subunit beta

Catalog No: #AP84431



Package Size: #AP84431-1 50ug #AP84431-2 100ug #AP84431-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-2 receptor subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IL-2 receptor subunit beta,High affinity IL-2 receptor subunit beta,p70-75,p75,CD122
Accession No.	Uniprot ID:P14784
Uniprot	P14784
GeneID	3560;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL2RB

Note: This product is for in vitro research use only