

Acetolactate synthase-like protein

Catalog No: #AP84432



Package Size: #AP84432-1 50ug #AP84432-2 100ug #AP84432-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Acetolactate synthase-like protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IlvB-like protein,AHAS
Accession No.	Uniprot ID:A1L0T0
Uniprot	A1L0T0
GeneID	10994;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ILVBL

Note: This product is for in vitro research use only