

Insulin-like peptide INSL6

Catalog No: #AP84433



Package Size: #AP84433-1 50ug #AP84433-2 100ug #AP84433-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Insulin-like peptide INSL6
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Insulin-like peptide 6,Relaxin/insulin-like factor 1,RIF1
Accession No.	Uniprot ID:Q9Y581
Uniprot	Q9Y581
GeneID	11172;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:INSL6

Note: This product is for in vitro research use only