

Insulin-like growth factor-binding protein 2

Catalog No: #AP84434



Package Size: #AP84434-1 50ug #AP84434-2 100ug #AP84434-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Insulin-like growth factor-binding protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	IBP-2,BRL-BP,Igfbp-2
Accession No.	Uniprot ID:P12843
Uniprot	P12843
GeneID	25662;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Igfbp2

Note: This product is for in vitro research use only