

Isocitrate dehydrogenase [NADP], mitochondrial

Catalog No: #AP84435



Package Size: #AP84435-1 50ug #AP84435-2 100ug #AP84435-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Isocitrate dehydrogenase [NADP], mitochondrial |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | IDH,ICD-M,IDP,NADP(+)-specific ICDH,Oxalosuccinate decarboxylase |
| Accession No. | Uniprot ID:P56574 |
| Uniprot | P56574 |
| GeneID | 361596; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ldh2

Note: This product is for in vitro research use only