

Sperm protamine P1

Catalog No: #AP84440



Package Size: #AP84440-1 50ug #AP84440-2 100ug #AP84440-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Sperm protamine P1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cysteine-rich protamine,Prm-1
Accession No.	Uniprot ID:P10118
Uniprot	P10118
GeneID	24685;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Prm1

Note: This product is for in vitro research use only