

Keratin, type I cytoskeletal 20

Catalog No: #AP84451

Package Size: #AP84451-1 50ug #AP84451-2 100ug #AP84451-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Keratin, type I cytoskeletal 20
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cytokeratin-20,CK-20,Keratin-20,K20,Protein IT
Accession No.	Uniprot ID:P35900
Uniprot	P35900
GeneID	54474;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:KRT20

Note: This product is for in vitro research use only