

Keratin, type I cytoskeletal 13

Catalog No: #AP84455



Package Size: #AP84455-1 50ug #AP84455-2 100ug #AP84455-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Keratin, type I cytoskeletal 13
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	47 kDa cytokeratin,Cytokeratin-13,CK-13,Keratin-13,K13,Krt1-13
Accession No.	Uniprot ID:P08730
Uniprot	P08730
GeneID	16663;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Krt13

Note: This product is for in vitro research use only