

Integrin alpha-5

Catalog No: #AP84459



Package Size: #AP84459-1 50ug #AP84459-2 100ug #AP84459-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Integrin alpha-5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CD49 antigen-like family member E,Fibronectin receptor subunit alpha,Integrin alpha-F,VLA-5,CD49e,FNRA
Accession No.	Uniprot ID:P08648
Uniprot	P08648
GeneID	3678;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ITGA5

Note: This product is for in vitro research use only