

Inosine triphosphate pyrophosphatase

Catalog No: #AP84463



Package Size: #AP84463-1 50ug #AP84463-2 100ug #AP84463-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Inosine triphosphate pyrophosphatase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	ITPase,Non-canonical purine NTP pyrophosphatase,Non-standard purine NTP pyrophosphatase,Nucleoside-triphosphate diphosphatase,Nucleoside-triphosphate pyrophosphatase,NTPase
Accession No.	Uniprot ID:F1NLH9
Uniprot	F1NLH9
GeneID	424390;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ITPA

Note: This product is for in vitro research use only