

## E3 ubiquitin-protein ligase HUWE1

Catalog No: #AP84465



Package Size: #AP84465-1 50ug #AP84465-2 100ug #AP84465-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	E3 ubiquitin-protein ligase HUWE1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	E3Histone,HECT, UBA and WWE domain-containing protein 1,HECT-type E3 ubiquitin transferase HUWE1,Upstream regulatory element-binding protein 1,URE-B1,Kiaa0312,Ureb1
Accession No.	Uniprot ID:Q7TMY8
Uniprot	Q7TMY8
GeneID	59026;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Huwe1

Note: This product is for in vitro research use only