

Integrin beta-1

Catalog No: #AP84467



Package Size: #AP84467-1 50ug #AP84467-2 100ug #AP84467-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Integrin beta-1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Fibronectin receptor subunit beta,VLA-4 subunit beta,CD29
Accession No.	Uniprot ID:P53712
Uniprot	P53712
GeneID	281876;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ITGB1

Note: This product is for in vitro research use only