

Tyrosine-protein kinase JAK3

Catalog No: #AP84470



Package Size: #AP84470-1 50ug #AP84470-2 100ug #AP84470-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tyrosine-protein kinase JAK3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Janus kinase 3, JAK-3, Leukocyte janus kinase, L-JAK
Accession No.	Uniprot ID:P52333
Uniprot	P52333
GeneID	3718;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: JAK3

Note: This product is for in vitro research use only