

Heat shock factor protein 4

Catalog No: #AP84471



Package Size: #AP84471-1 50ug #AP84471-2 100ug #AP84471-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Heat shock factor protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	HSF 4,Heat shock transcription factor 4,HSTF 4
Accession No.	Uniprot ID:Q9R0L1
Uniprot	Q9R0L1
GeneID	26386;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Hsf4

Note: This product is for in vitro research use only