

Interferon-induced protein with tetratricopeptide repeats 1

Catalog No: #AP84477

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP84477-1 50ug #AP84477-2 100ug #AP84477-3 1mg

Description

Product Name	Interferon-induced protein with tetratricopeptide repeats 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IFIT-1,Glucocorticoid-attenuated response gene 16 protein,GARG-16,Interferon-induced 56 kDa protein,IFI-56K,Garg16,lfi56,lsg56
Accession No.	Uniprot ID:Q64282
Uniprot	Q64282
GeneID	15957;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:lfit1

Note: This product is for in vitro research use only