

Insulin-like growth factor II

Catalog No: #AP84479



Package Size: #AP84479-1 50ug #AP84479-2 100ug #AP84479-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Insulin-like growth factor II
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IGF-II,Somatomedin-A,T3M-11-derived growth factor,PP1446
Accession No.	Uniprot ID:P01344
Uniprot	P01344
GeneID	3481;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IGF2

Note: This product is for in vitro research use only