

Type I iodothyronine deiodinase

Catalog No: #AP84482



Package Size: #AP84482-1 50ug #AP84482-2 100ug #AP84482-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Type I iodothyronine deiodinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	5DI,DIOI,Type 1 DI,Type-I 5'-deiodinase,ITDI1,TXDI1
Accession No.	Uniprot ID:P49894
Uniprot	P49894
GeneID	403635;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:DIO1

Note: This product is for in vitro research use only