

Ras GTPase-activating-like protein IQGAP1

Catalog No: #AP84484



Package Size: #AP84484-1 50ug #AP84484-2 100ug #AP84484-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ras GTPase-activating-like protein IQGAP1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Accession No.	Uniprot ID:Q9JKF1
Uniprot	Q9JKF1
GeneID	29875;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Iqgap1

Note: This product is for in vitro research use only