

Interferon alpha/beta receptor 1

Catalog No: #AP84493



Package Size: #AP84493-1 50ug #AP84493-2 100ug #AP84493-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interferon alpha/beta receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IFN-R-1,Cytokine receptor class-II member 1,Cytokine receptor family 2 member 1,CRF2-1,Type I interferon receptor 1,IFNAR
Accession No.	Uniprot ID:P17181
Uniprot	P17181
GeneID	3454;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IFNAR1

Note: This product is for in vitro research use only