

Inorganic pyrophosphatase

Catalog No: #AP84495



Package Size: #AP84495-1 50ug #AP84495-2 100ug #AP84495-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Inorganic pyrophosphatase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Pyrophosphate phospho-hydrolase,PPase,Pp,Pyp
Accession No.	Uniprot ID:Q9D819
Uniprot	Q9D819
GeneID	67895;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ppa1

Note: This product is for in vitro research use only