

Ras GTPase-activating-like protein IQGAP1

Catalog No: #AP84496



Package Size: #AP84496-1 50ug #AP84496-2 100ug #AP84496-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ras GTPase-activating-like protein IQGAP1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	p195,KIAA0051
Accession No.	Uniprot ID:P46940
Uniprot	P46940
GeneID	8826;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IQGAP1

Note: This product is for in vitro research use only