

Lysosome-associated membrane glycoprotein 2

Catalog No: #AP84500



Package Size: #AP84500-1 50ug #AP84500-2 100ug #AP84500-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lysosome-associated membrane glycoprotein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	LAMP-2,CD107 antigen-like family member B,LGP-110,LGP-96,Lysosomal membrane glycoprotein type B,LGP-B,CD107b,Lamp-2
Accession No.	Uniprot ID:P17046
Uniprot	P17046
GeneID	24944;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lamp2

Note: This product is for in vitro research use only