

Lysosomal-associated transmembrane protein 4B

Catalog No: #AP84501



Package Size: #AP84501-1 50ug #AP84501-2 100ug #AP84501-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lysosomal-associated transmembrane protein 4B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Lysosome-associated transmembrane protein 4-beta,PSEC0001
Accession No.	Uniprot ID:Q86V14
Uniprot	Q86V14
GeneID	55353;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: LAPT4B

Note: This product is for in vitro research use only