

Inhibitor of nuclear factor kappa-B kinase subunit epsilon

Catalog No: #AP84504

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP84504-1 50ug #AP84504-2 100ug #AP84504-3 1mg

Description

Product Name	Inhibitor of nuclear factor kappa-B kinase subunit epsilon
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	I-kappa-B kinase epsilon,Inducible I kappa-B kinase,IKK-i,IKKE,IKKI,KIAA0151
Accession No.	Uniprot ID:Q14164
Uniprot	Q14164
GeneID	9641;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IKBKE

Note: This product is for in vitro research use only