

Interleukin-2 receptor subunit alpha

Catalog No: #AP84505



Package Size: #AP84505-1 50ug #AP84505-2 100ug #AP84505-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Interleukin-2 receptor subunit alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IL-2 receptor subunit alpha,TAC antigen,p55,CD25
Accession No.	Uniprot ID:P01589
Uniprot	P01589
GeneID	3559;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL2RA

Note: This product is for in vitro research use only