

Cytokine receptor common subunit gamma

Catalog No: #AP84506



Package Size: #AP84506-1 50ug #AP84506-2 100ug #AP84506-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cytokine receptor common subunit gamma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Interleukin-2 receptor subunit gamma,IL-2 receptor subunit gamma,gammaC,p64,CD132
Accession No.	Uniprot ID:Q95118
Uniprot	Q95118
GeneID	281862;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IL2RG

Note: This product is for in vitro research use only