

Interferon regulatory factor 2-binding protein 2

Catalog No: #AP84507



Package Size: #AP84507-1 50ug #AP84507-2 100ug #AP84507-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interferon regulatory factor 2-binding protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	IRF-2-binding protein 2
Accession No.	Uniprot ID:Q7Z5L9
Uniprot	Q7Z5L9
GeneID	359948;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IRF2BP2

Note: This product is for in vitro research use only