

# 6-phosphofructokinase, muscle type

Catalog No: #AP84509



Package Size: #AP84509-1 50ug #AP84509-2 100ug #AP84509-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	6-phosphofructokinase, muscle type
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Phosphofructo-1-kinase isozyme A,PFK-A,Phosphofructokinase-M,Phosphofructokinase 1,Phosphohexokinase,PFKX
Accession No.	Uniprot ID:P08237
Uniprot	P08237
GeneID	5213;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:PFKM

Note: This product is for in vitro research use only