

UMP-CMP kinase

Catalog No: #AP84511

Package Size: #AP84511-1 50ug #AP84511-2 100ug #AP84511-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	UMP-CMP kinase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Deoxycytidylate kinase,CK,Nucleoside-diphosphate kinase,Uridine monophosphate/cytidine monophosphate kinase,UMP/CMP kinase,Cmk,Cmpk,Uck,Umk,Umpk
Accession No.	Uniprot ID:Q9DBP5
Uniprot	Q9DBP5
GeneID	66588;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cmpk1

Note: This product is for in vitro research use only