

## Keratin, type I cytoskeletal 17

Catalog No: #AP84518



Package Size: #AP84518-1 50ug #AP84518-2 100ug #AP84518-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Keratin, type I cytoskeletal 17
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	39.1,Cytokeratin-17,CK-17,Keratin-17,K17
Accession No.	Uniprot ID:Q04695
Uniprot	Q04695
GeneID	3872;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:KRT17

Note: This product is for in vitro research use only