

## Keratin, type II cytoskeletal 6A

Catalog No: #AP84520



Package Size: #AP84520-1 50ug #AP84520-2 100ug #AP84520-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Keratin, type II cytoskeletal 6A
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Cytokeratin-6A,CK-6A,Keratin-6-alpha,mK6-alpha,Keratin-6A,K6A,Ker2,Krt2-6,Krt2-6a,Krt6
Accession No.	Uniprot ID:P50446
Uniprot	P50446
GeneID	16687;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Krt6a

Note: This product is for in vitro research use only