

Hyaluronan synthase 2

Catalog No: #AP84521



Package Size: #AP84521-1 50ug #AP84521-2 100ug #AP84521-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Hyaluronan synthase 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Hyaluronate synthase 2, Hyaluronic acid synthase 2, HA synthase 2
Accession No.	Uniprot ID:O97711
Uniprot	O97711
GeneID	281220;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: HAS2

Note: This product is for in vitro research use only