

Insulin-like growth factor-binding protein 4

Catalog No: #AP84522



Package Size: #AP84522-1 50ug #AP84522-2 100ug #AP84522-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Insulin-like growth factor-binding protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	IBP-4
Accession No.	Uniprot ID:Q05716
Uniprot	Q05716
GeneID	282262;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: IGFBP4

Note: This product is for in vitro research use only