

Isocitrate dehydrogenase [NADP], mitochondrial

Catalog No: #AP84523



Package Size: #AP84523-1 50ug #AP84523-2 100ug #AP84523-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Isocitrate dehydrogenase [NADP], mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IDH,ICD-M,IDP,NADP(+)-specific ICDH,Oxalosuccinate decarboxylase
Accession No.	Uniprot ID:P54071
Uniprot	P54071
GeneID	269951;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ldh2

Note: This product is for in vitro research use only