

Keratin, type I cytoskeletal 10

Catalog No: #AP84525



Package Size: #AP84525-1 50ug #AP84525-2 100ug #AP84525-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Keratin, type I cytoskeletal 10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cytokeratin-10,CK-10,Keratin-10,K10,KPP
Accession No.	Uniprot ID:P13645
Uniprot	P13645
GeneID	3858;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:KRT10

Note: This product is for in vitro research use only