

# Keratin, type I cytoskeletal 19

Catalog No: #AP84526



Package Size: #AP84526-1 50ug #AP84526-2 100ug #AP84526-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Keratin, type I cytoskeletal 19
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cytokeratin-19,CK-19,Keratin-19,K19,Type I keratin Ka19,Ka19,Krt1-19
Accession No.	Uniprot ID:Q63279
Uniprot	Q63279
GeneID	360626;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Krt19

Note: This product is for in vitro research use only