

Golgi apparatus protein 1

Catalog No: #AP84529



Package Size: #AP84529-1 50ug #AP84529-2 100ug #AP84529-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Golgi apparatus protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	Cysteine-rich fibroblast growth factor receptor,CFR
Accession No.	Uniprot ID:Q02391
Uniprot	Q02391
GeneID	396492;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:GLG1

Note: This product is for in vitro research use only